

SEAMLESS INTEGRATION

Import images directly into patient files

Magewell Image Capture Kit

A simple solution for creating digital stills from a live video stream. Easily archive files from your microscope's camera to a computer during an exam.



- One-click image or video capture to your PC using Capture Express software.
- No need to transfer files using an external media card.
- Improve workflow by real-time capture to patient software that recognizes an external camera or webcam.
- Simultaneously stream live video to your computer and multiple monitors within your operatory.
- Versatile Air Mouse remote control allows you to capture images untethered from your computer.
- Use with most existing equipment; compatible with newer monitors, touch screen displays, and computers that meet minimum requirements.
- Plug and play – no drivers required to transfer files.
- Minimal configuration for connecting to patient software.*



*Due to the wide variety of programs, additional support from your IT group or from the third-party software provider may occasionally be needed.

GLOBAL
SURGICAL™ CORPORATION

Magewell Image Capture Kit



The Magewell Kit also contains an owner's manual to assist with setup and configuration variables. In it, you will find helpful tips for installing and using the software.

MINIMUM REQUIREMENTS:

- 1080P HD camera with HDMI cable
- Windows 7, 8, or 10 OS
- i5 Processor
- USB 3.0 port

COMPATIBLE:

- 1080P HD monitor
- Patient software that recognizes external camera or webcam through USB Video Class (UVC) driver.

INCLUDED COMPONENTS:

- Magewell Capture Device:
- 2' USB 3.0 cable
- Capture Express Software on USB flash drive
- Wireless Air Mouse
- Owner's Manual