

# GLOBAL QUALITY. NOW IN 4K.

## HANDS-FREE 4K VIDEO & HD PHOTO SYSTEM

Supercharge your documentation, case presentation and patient education with crystal-clear 4K video and HD photography. Connects easily to your Global microscope and includes a foot pedal for hands-free operation, so you can stay focused - with no interruptions to your treatment.

- Quick, on-the-fly video and image capture
- Foot pedal for hands-free operation
- Easily connects via a C-mount adapter
- USB integration
- 4K/1080p/720p video @ 30 fps
- 8.3 MP images

